



Anti-SNAPC1 (internal region) polyclonal antibody (DPABH-27248)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Part of the SNAPc complex required for the transcription of both RNA polymerase II and III small-nuclear RNA genes. Binds to the proximal sequence element (PSE), a non-TATA-box basal promoter element common to these 2 types of genes. Recruits TBP and BRF2 to the U6 snRNA TATA box.
Immunogen	Synthetic peptide, derived from an internal region of Mouse SNAP43 (UniProt ID: Q8K0S9).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse
Purification	Whole antiserum
Conjugate	Unconjugated
Applications	WB, ICC
Format	Liquid
Size	200 µl
Buffer	Constituent: Whole serum
Preservative	None
Storage	Store at 4°C short term (1-2 weeks). Aliquot and store at -20°C long term. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Gene Name	SNAPC1 small nuclear RNA activating complex, polypeptide 2 [Mus musculus]
Official Symbol	SNAPC1
Synonyms	SNAPC1; small nuclear RNA activating complex, polypeptide 1; AU015787; 2700033G17Rik; snRNA-activating protein complex subunit 1; SNAPc subunit 1; SNAPc 43 kDa subunit; snRNA-activating protein complex 43 kDa subunit; small nuclear RNA-activating complex polypeptide 1;
Entrez Gene ID	75627
Protein Refseq	NP_848479.1
UniProt ID	Q8K0S9
Pathway	Gene Expression; RNA Polymerase III Abortive And Retractive Initiation; RNA Polymerase III Transcription Initiation;
Function	DNA binding; molecular_function;