



Anti-MAP3K3 (N-terminal) polyclonal antibody (DPABH-27255)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Component of a protein kinase signal transduction cascade. Mediates activation of the NF-kappa-B, AP1 and DDIT3 transcriptional regulators.
Immunogen	Synthetic peptide derived from sequence near the N-terminus of Mouse MEKK3.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse
Purification	Whole antiserum
Conjugate	Unconjugated
Applications	WB, IHC-P
Format	Liquid
Size	50 µl
Buffer	Constituent: Whole serum
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name [MAP3K3 mitogen-activated protein kinase kinase kinase 4 \[Mus musculus \]](#)

Official Symbol	MAP3K3
Synonyms	MAP3K3; mitogen-activated protein kinase kinase kinase 3; Mekk3; MAPKKK3; AW548911; mKIAA4031; MEKK 3; MEK kinase 3; MAPK/ERK kinase kinase 3; mitogen activated protein kinase kinase kinase 3;
Entrez Gene ID	26406
Protein Refseq	NP_036077.1
UniProt ID	Q61084
Pathway	Cytokine Signaling in Immune system; EGFR1 Signaling Pathway; GnRH signaling pathway; HTLV-I infection
Function	ATP binding; MAP kinase kinase kinase activity; kinase activity; metal ion binding