



Anti-TRIM2 (aa 233-247) polyclonal antibody (DPABH-27348)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	May contribute to the alteration of neural cellular mechanisms.
Specificity	DPABH-27348 is expected to recognize both reported isoforms (NP_056086.2; NP_001123539.1).
Immunogen	Synthetic peptide: C-STFDELQKTLNVRKS, corresponding to internal sequence amino acids 233-247 of Mouse RNF86 (NP_056086.2; NP_001123539.1).
Isotype	IgG
Source/Host	Goat
Species Reactivity	Mouse
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Size	200 µl
Buffer	pH: 7.30; Constituents: 99% Tris buffered saline, 0.5% BSA
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Gene Name	TRIM2 tripartite motif-containing 3 [Mus musculus]
Official Symbol	TRIM2
Synonyms	TRIM2; tripartite motif-containing 2; narf; mKIAA0517; tripartite motif-containing protein 2; tripartite motif protein 2; tripartite motif protein TRIM2; E3 ubiquitin-protein ligase TRIM2; neural activity-related ring finger protein;
Entrez Gene ID	80890
Protein Refseq	NP_001258654.1
UniProt ID	E9QKC6
Function	UDP-N-acetylmuramoylalanyl-D-glutamyl-2,6-diaminopimelate-D-alanyl-D-alanine ligase activity; coenzyme F420-0 gamma-glutamyl ligase activity; coenzyme F420-2 alpha-glutamyl ligase activity; ligase activity