



Anti-RHOX8 (aa 201-250) polyclonal antibody (DPABH-27586)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	The function of this protein remains unknown.
Immunogen	Synthetic peptide corresponding to a region within internal sequence amino acids 201-250 (EAAAAAARDETTAGSSAPVDDRSHDAGASSNGEDRGQGEELIPGGTKG L) of Mouse Rhox8 (NP_001004193).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Size	50 µg
Buffer	Constituents: 2% Sucrose, PBS
Preservative	None
Storage	Store at 4°C short term (1-2 weeks). Aliquot and store at -20°C long term. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Gene Name	RHOX8 reproductive homeobox 9 [Mus musculus]
Official Symbol	RHOX8
Synonyms	RHOX8; reproductive homeobox 8; Tox; reproductive homeobox on X chromosome, 8; homeoprotein TOX;
Entrez Gene ID	434768
Protein Refseq	NP_001004193.1
UniProt ID	Q6VSS7
Function	DNA binding;