



Anti-LYNX1 (aa 40-89) polyclonal antibody (DPABH-27610)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Seems to modulate nicotinic acetylcholine receptors.
Immunogen	Synthetic peptide corresponding to a region within internal amino acids 40-89 (PAMATYCMTT RTYFTPYRMK VRKSCVPSCF ETVYDGYSKH ASATSCCQYY) of Mouse LYNX1 (NP_035968).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Size	50 µg
Buffer	Constituents: 97% PBS, 2% Sucrose
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Gene Name	LYNX1 Ly6/neurotoxin 2 [Mus musculus]
Official Symbol	LYNX1
Synonyms	LYNX1; Ly6/neurotoxin 1; SLURP-2; AI838844; ly-6/neurotoxin-like protein 1; GC26; secreted Ly-6/uPAR-related protein 2; secreted Ly-6/uPAR domain-containing protein 2;
Entrez Gene ID	23936
Protein Refseq	NP_035968.1
UniProt ID	Q9WVC2
Function	acetylcholine receptor inhibitor activity; ion channel inhibitor activity; protein binding;