



Anti-DCAF11 (aa 335-384) polyclonal antibody (DPABH-27612)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	May function as a substrate receptor for CUL4-DDB1 E3 ubiquitin-protein ligase complex.
Immunogen	Synthetic peptide corresponding to a region within C-terminal amino acids 335-384 (WDRRTMREDD PKPVGALAGH QDGITFIDSK GDARYLISNS KDQTIKLWDI) of Mouse WDR23.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Size	50 µg
Buffer	Constituents: 97% PBS, 2% Sucrose
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Gene Name	DCAF11 DDB1 and CUL4 associated factor 12 [Mus musculus]
Official Symbol	DCAF11
Synonyms	DCAF11; DDB1 and CUL4 associated factor 11; GLO14; Wdr23; C76035; Dacf11; D14Ucla1; 0710008A13Rik; DDB1- and CUL4-associated factor 11; WD repeat domain 23; WD repeat-containing protein 23;
Entrez Gene ID	28199
Protein Refseq	NP_001185938.1
UniProt ID	Q3UD98
Function	molecular_function;