



Rabbit anti-Mouse Cathepsin Z polyclonal antibody (DPABH-27623)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Exhibits carboxy-monopeptidase as well as carboxy-dipeptidase activity.
Immunogen	Synthetic peptide corresponding to a region within N terminal amino acids 21-70 (SAARARLYFR SGQTCYHPIR GDQLALLGRR TYPRPHEYLS PADLPKNWDW) of Mouse Cathepsin Z.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Size	50 µg
Buffer	Constituents: 97% PBS, 2% Sucrose
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Gene Name	CTSZ cathepsin Z [Mus musculus]
Official Symbol	CTSZ
Synonyms	CTSZ; cathepsin Z; CTSX; AI787083; AU019819; D2Wsu143e; cathepsin X;
Entrez Gene ID	64138
Protein Refseq	NP_071720.1
UniProt ID	Q545I6
Pathway	Clathrin derived vesicle budding; Lysosome; Membrane Trafficking; Metabolism of proteins
Function	cysteine-type peptidase activity; hydrolase activity; peptidase activity;