



Anti-MUC19 (aa 7461-7474) polyclonal antibody (DPABH-27826)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	May function in ocular mucus homeostasis.
Immunogen	Synthetic peptide: ECKRSVKYNYETFQ, corresponding to C terminal amino acids 7461-7474 of Mouse Muc19 (NP_997126.2).
Isotype	IgG
Source/Host	Goat
Species Reactivity	Mouse
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	ICC/IF, ELISA
Format	Liquid
Size	200 µl
Buffer	pH: 7.30; Constituents: 99% Tris buffered saline, 0.5% BSA
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Gene Name [MUC19 mucin 20 \[Mus musculus \]](#)

Official Symbol	MUC19
Synonyms	MUC19; mucin 19; sld; MUC-19; mucin-19; sublingual apomucin; gel-forming secreted mucin-19;
Entrez Gene ID	239611
Protein Refseq	NP_997126.2
UniProt ID	Q6PZE0
Pathway	Metabolism of proteins; O-linked glycosylation of mucins; Termination of O-glycan biosynthesis;
Function	molecular_function;