



Anti-BPIFA2 (aa 146-160) polyclonal antibody (DPABH-27830)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	BPIFA2 (BPI fold containing family A, member 2) is a protein-coding gene. GO annotations related to this gene include lipid binding.
Immunogen	Synthetic peptide: KTNAQTGLPEVTIGK, corresponding to internal sequence amino acids 146-160 of Mouse C20orf70.
Isotype	IgG
Source/Host	Goat
Species Reactivity	Mouse
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	ICC/IF, ELISA
Format	Liquid
Size	200 µl
Buffer	pH: 7.30; Constituents: 99.48% Tris buffered saline, 0.5% BSA
Preservative	0.02% Sodium Azide
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name	BPIFA2 BPI fold containing family A, member 3 [Mus musculus]
Official Symbol	BPIFA2
Synonyms	BPIFA2; BPI fold containing family A, member 2; Psp; Bpifa2e; BPI fold-containing family A member 2; parotid secretory protein; BPI fold containing family A, member 2E;
Entrez Gene ID	19194
Protein Refseq	NP_032979.1
UniProt ID	P07743
Function	lipid binding; lipopolysaccharide binding;