



# Anti-FAM49B (aa 2-134) polyclonal antibody (DPABH-28004)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Antigen Description</b>	FAM49B (family with sequence similarity 49, member B) is a protein-coding gene. Diseases associated with FAM49B include lung cancer. An important paralog of this gene is FAM49A.
<b>Immunogen</b>	Recombinant fragment, corresponding to amino acids 2-134 of Human FAM49B (Q9NUQ9).
<b>Isotype</b>	IgG
<b>Source/Host</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Purification</b>	Immunogen affinity purified
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	ICC/IF, IHC-P, WB
<b>Format</b>	Liquid
<b>Size</b>	100 µl
<b>Buffer</b>	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
<b>Preservative</b>	0.02% Sodium Azide
<b>Storage</b>	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

## GENE INFORMATION

**Gene Name** [FAM49B family with sequence similarity 50, member B \[ Mus musculus \]](#)

<b>Official Symbol</b>	FAM49B
<b>Synonyms</b>	FAM49B; family with sequence similarity 49, member B; AW122079; 0910001A06Rik; protein FAM49B;
<b>Entrez Gene ID</b>	<a href="#">223601</a>
<b>Protein Refseq</b>	<a href="#">NP_659095.1</a>
<b>UniProt ID</b>	<a href="#">Q921M7</a>
<b>Function</b>	protein binding;