



Anti-DYTN (internal region) polyclonal antibody (DPABH-28212)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	DYTN (dystrotelin) is a protein-coding gene. GO annotations related to this gene include zinc ion binding and calcium ion binding. An important paralog of this gene is DTNA.
Immunogen	Synthetic peptide derived from internal domain of Mouse Dystrotelin.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse
Purification	Whole antiserum
Conjugate	Unconjugated
Applications	WB, ELISA, ICC
Format	Liquid
Size	200 µl
Buffer	Constituents: Whole serum
Preservative	None
Storage	Store at 4°C short term (1-2 weeks). Aliquot and store at -20°C long term. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Gene Name	DYTN dystrotelin [Mus musculus]
Official Symbol	DYTN
Synonyms	DYTN; dystrotelin; Gm215;
Entrez Gene ID	241073
Protein Refseq	NP_001075127.1
UniProt ID	A2CI98
Function	calcium ion binding; metal ion binding; molecular_function; zinc ion binding