



Anti-TLX1 (aa 1-100) polyclonal antibody (DPABH-28273)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Controls the genesis of the spleen. Binds to the DNA sequence 5-GGCGGTAAGTGG-3.
Immunogen	Synthetic peptide conjugated to KLH derived from within residues 1 - 100 of Mouse HOX11.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Size	100 µg
Buffer	pH: 7.40; Constituent: PBS
Preservative	0.02% Sodium Azide
Storage	Store at 4°C short term (1-2 weeks). Aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Gene Name [TLX1 T cell leukemia, homeobox 2 \[Mus musculus \]](#)

Official Symbol	TLX1
Synonyms	TLX1; T cell leukemia, homeobox 1; Hox11; Hox-11; T-cell leukemia homeobox protein 1; homeobox TLX-1; homeobox protein Hox-11; T-cell leukemia, homeobox 1;
Entrez Gene ID	21908
Protein Refseq	NP_068701.1
UniProt ID	Q9QX99
Pathway	Transcriptional misregulation in cancer;
Function	DNA binding; protein heterodimerization activity;