



Anti-BOLA2 (aa 1-50) polyclonal antibody (DPABH-28283)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	BOLA2 (boIA family member 2) is a protein-coding gene. An important paralog of this gene is BOLA2B.
Immunogen	Synthetic peptide corresponding to a region within N terminal amino acids 1-50 (MELSADYLRE KLRQDLEAEH VEVEDTTLNR CATSFRVLVV SAKFEGKPLL) of Mouse BOLA2 (NP_780312)
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Size	50 µg
Buffer	Constituents: 2% Sucrose, PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Gene Name	BOLA2 bolA-like 3 (E. coli) [Mus musculus]
Official Symbol	BOLA2
Synonyms	BOLA2; bolA-like 2 (E. coli); BolA-2; 1110025L05Rik; bolA-like protein 2;
Entrez Gene ID	66162
Protein Refseq	NP_780312.1
UniProt ID	Q8BGS2
Function	molecular_function;