



Anti-NKX2-4 (aa 301-350) polyclonal antibody (DPABH-28372)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Probable transcription factor.
Immunogen	Synthetic peptide corresponding to a region within the C terminal amino acids 301-350 (QQPQAPTPAP ELEELSPSP ALHGPGGGLA ALDAATGDYG GGVLGANLLY) of Mouse Nkx2.4 (NP_075993).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Size	50 µg
Buffer	Constituents: 2% Sucrose, PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Gene Name	NKX2-4 NK2 homeobox 5 [Mus musculus]
Official Symbol	NKX2-4
Synonyms	NKX2-4; NK2 homeobox 4; tinman; Nkx-2.4; 1700001P03Rik; homeobox protein Nkx-2.4; homeobox protein NK-2 homolog D; NK2 transcription factor related, locus 4;
Entrez Gene ID	228731
Protein Refseq	NP_075993.1
UniProt ID	Q9EQM3
Function	DNA binding; sequence-specific DNA binding; sequence-specific DNA binding transcription factor activity;