



Anti-PCID2 (aa 169-218) polyclonal antibody (DPABH-28815)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	PCID2 (PCI domain containing 2) is a protein-coding gene.
Immunogen	Synthetic peptide corresponding to a region within C terminal amino acids 169-218 (MLFLVNQLFK IYFKINKLHL CKPLIRAIDS SNLKDDYSTA QRVTYKYYVG) of Rat RGD1307041 (NP_001162615).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Rat
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Size	50 µg
Buffer	Constituents: PBS, 2% Sucrose
Preservative	None
Storage	Store at 4°C. Aliquot and store at -20°C long-term.

GENE INFORMATION

Gene Name	PCID2 PCI domain containing 3 [Rattus norvegicus]
Official Symbol	PCID2
Synonyms	PCID2; PCI domain containing 2; RGD1307041; uncharacterized protein LOC361182;
Entrez Gene ID	361182
Protein Refseq	NP_001162615.1
UniProt ID	B1H263