



Anti-FRMD5 (aa 164-213) polyclonal antibody (DPABH-28920)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	FRMD5 (FERM domain containing 5) is a protein-coding gene. GO annotations related to this gene include cytoskeletal protein binding and molecular_function. An important paralog of this gene is EPB41L5.
Immunogen	Synthetic peptide corresponding to a region within N terminal amino acids 164-213 (FPKHSEKLEK KIAEIHKTEL SGQTPATSEL NFLRKAQTLE TYGVDPHPCK) of Rat FRMD5 (XP_230513).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Rat
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Size	50 µg
Buffer	Constituents: 97% PBS, 2% Sucrose
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Gene Name	FRMD5 FERM domain containing 6 [Rattus norvegicus]
Official Symbol	FRMD5
Synonyms	FRMD5; FERM domain containing 5; RGD1561554;
Entrez Gene ID	311362
UniProt ID	F1LX12
Function	cytoskeletal protein binding;