



Anti-PYY polyclonal antibody (DPABH-29210)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	This gut peptide inhibits exocrine pancreatic secretion, has a vasoconstrictory action and inhibits jejunal and colonic mobility.
Immunogen	Synthetic peptide conjugated to KLH derived from within residues 1 to the C-terminus of Rat Peptide YY.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse, Rat, Human, Marmoset
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	IHC-P
Format	Liquid
Size	100 µg
Buffer	pH: 7.40; Constituent: PBS
Preservative	0.02% Sodium Azide
Storage	Store at 4°C short term (1-2 weeks). Aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Gene Name [PYY peptide YY \(mapped\) \[Rattus norvegicus \]](#)

Official Symbol	PYY
Synonyms	PYY; peptide YY (mapped); Yy; GHYY; RATGHYY; peptide-YY; peptide YY; peptide tyrosine tyrosine; peptide tyrosine-tyrosine (YY);
Entrez Gene ID	287730
Protein Refseq	NP_001029252.1
UniProt ID	P10631
Function	neuropeptide hormone activity; neuropeptide hormone activity; neuropeptide hormone activity;