



Anti-MAEL polyclonal antibody (DPABT-H11346)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Anti-MAEL Polyclonal Antibody
Target	MAEL
Immunogen	Recombinant protein corresponding to amino acids of human MAEL.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Antigen affinity purification
Conjugate	Unconjugated
Applications	WB, IHC
Sequence Similarities	PVFTPLRRPGMLVPKQNVSPDMSALSLKGDQALLGGIFYFLNIFSHGELPPHCEQRFLP CEIGCVKYSLQEGIMADFH
Format	Liquid
Buffer	In PBS, pH 7.5 (40% glycerol, 0.02% sodium azide)
Preservative	0.02% Sodium Azide
Storage	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name	MAEL maelstrom homolog (Drosophila) [Homo sapiens]
Official Symbol	MAEL
Synonyms	MAEL; maelstrom homolog (Drosophila); protein maelstrom homolog; cancer/testis antigen 128; CT128; FLJ14904; RP11-102C16.1;
Entrez Gene ID	84944
Protein Refseq	NP_116247
UniProt ID	Q96JY0
Chromosome Location	1q24.1
Function	DNA binding;