



Anti-PKHD1 polyclonal antibody (DPABT-H12329)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Anti-PKHD1 Polyclonal Antibody
Target	PKHD1
Immunogen	Recombinant protein corresponding to amino acids of human PKHD1.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Antigen affinity purification
Conjugate	Unconjugated
Applications	IHC
Sequence Similarities	KTTTVNYVRDTLSNPRGWMALLLDQETYSLQSENLWINRSLQYSATFDNFAPGNLYLLLVH TDLPPYPDILLRCGSRVGLSFPFLSPGQNNQ
Format	Liquid
Size	100 µl
Buffer	In PBS, pH 7.5 (40% glycerol, 0.02% sodium azide)
Preservative	0.02% Sodium Azide
Storage	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name	PKHD1 polycystic kidney and hepatic disease2(autosomal recessive) [Homo sapiens]
Official Symbol	PKHD1
Synonyms	PKHD1; polycystic kidney and hepatic disease2(autosomal recessive); TIG multiple domains2, TIGM1; fibrocystin; ARPKD; FCYT; polyductin; tigmin; TIG multiple domains 1; polycystic kidney and hepatic disease2protein; TIGM1; FLJ46150; DKFZp686C01112;
Entrez Gene ID	5314
Protein Refseq	NP_619639
UniProt ID	P08F94
Chromosome Location	6p21.2-p12
Function	protein binding; receptor activity;