



Anti-PLXNA4 polyclonal antibody (DPABT-H12377)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Anti-PLXNA4 Polyclonal Antibody
Target	PLXNA4
Immunogen	Recombinant protein corresponding to amino acids of human PLXNA4.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Antigen affinity purification
Conjugate	Unconjugated
Applications	IHC
Sequence Similarities	SHVKVAGVECSPLVDGYIPAEQIVCEMGEAKPSQHAGFVEICVAVCRPEFMARSSQLYYF MTLTLSDLKPSRGPMSGGTQVTI
Format	Liquid
Buffer	In PBS, pH 7.5 (40% glycerol, 0.02% sodium azide)
Preservative	0.02% Sodium Azide
Storage	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name	PLXNA4 plexin A4 [Homo sapiens]
Official Symbol	PLXNA4
Synonyms	PLXNA4; plexin A4; plexin A4, A , plexin A4, B , PLXNA4A, PLXNA4B; DKFZp434G0625PRO34003; FAYV2820; KIAA1550; PLEXA4;
Entrez Gene ID	23655
Chromosome Location	7q32.3