



Anti-SLC15A3 polyclonal antibody (DPABT-H13167)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Anti-SLC15A3 Polyclonal Antibody
Target	SLC15A3
Immunogen	Recombinant protein corresponding to amino acids of human SLC15A3.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Antigen affinity purification
Conjugate	Unconjugated
Applications	IHC
Sequence Similarities	PMGSQVSSMLKLALQNCCPQLWQRHSARDRQCARVLADERSPQPGASPQEDIANFQ
Format	Liquid
Buffer	In PBS, pH 7.5 (40% glycerol, 0.02% sodium azide)
Preservative	0.02% Sodium Azide
Storage	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name	SLC15A3 solute carrier family 15, member4[Homo sapiens]
Official Symbol	SLC15A3
Synonyms	SLC15A3; solute carrier family 15, member 3; solute carrier family 15 member 3; hPTR3; PHT2; peptide transporter 3; osteoclast transporter; peptide/histidine transporter 2; OCTP; PTR3; FLJ26631;
Entrez Gene ID	51296
Protein Refseq	NP_057666
UniProt ID	Q8IY34
Chromosome Location	11q12.2
Pathway	Proton/oligonucleotide cotransporters, organism-specific biosystem; SLC-mediated transmembrane transport, organism-specific biosystem; Transmembrane transport of small molecules, organism-specific biosystem; Transport of inorganic cations/anions and amino acids/oligopeptides, organism-specific biosystem;
Function	symporter activity; transporter activity;