



Anti-GPR55 polyclonal antibody (DPABT-H17111)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Anti-GPR55 Polyclonal Antibody
Target	GPR55
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	N/A
Conjugate	Unconjugated
Applications	IHC
Preservative	None
Storage	-70 °C
Ship	4 °C

GENE INFORMATION

Gene Name	GPR55 G protein-coupled receptor 55 [Homo sapiens]
Official Symbol	GPR55
Synonyms	GPR55; G protein-coupled receptor 55; G-protein coupled receptor 55;
Entrez Gene ID	9290

Protein Refseq	NP_005674
UniProt ID	A8K858
Chromosome Location	2q37
Pathway	Class A/1 (Rhodopsin-like receptors), organism-specific biosystem; G alpha (i) signalling events, organism-specific biosystem; GPCR downstream signaling, organism-specific biosystem; GPCR ligand binding, organism-specific biosystem; GPCRs, Other, organism
Function	G-protein coupled receptor activity; cannabinoid receptor activity; receptor activity; signal transducer activity;