



Anti-PRKCD polyclonal antibody (DPABT-H20217)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Anti-PRKCD Polyclonal Antibody
Specificity	Recognizes the ~77 kDa PKCdelta protein.
Target	PRKCD
Immunogen	a synthetic peptide
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	N/A
Conjugate	Unconjugated
Applications	WB
Format	Purified
Size	500 µl
Preservative	None
Storage	Following initial thaw, aliquot and freeze (-70 °C). Avoid repeated freeze-thaw.

GENE INFORMATION

Gene Name [PRKCD protein kinase C, delta \[Homo sapiens \]](#)

Official Symbol	PRKCD
Synonyms	PRKCD; protein kinase C, delta; protein kinase C delta type; protein kinase C delta VIII; tyrosine-protein kinase PRKCD; MAY1; PKCD; nPKC-delta; MGC49908;
Entrez Gene ID	5580
Protein Refseq	NP_006245
UniProt ID	A0A024R328
Chromosome Location	3p21.31
Pathway	Alpha-synuclein signaling, organism-specific biosystem; Alpha6-Beta4 Integrin Signaling Pathway, organism-specific biosystem; Apoptosis, organism-specific biosystem; Apoptotic cleavage of cellular proteins, organism-specific biosystem; Apoptotic executionphase, organism-specific biosystem; B Cell Receptor Signaling Pathway, organism-specific biosystem; Ca-dependent events, organism-specific biosystem;
Function	ATP binding; calcium-independent protein kinase C activity; enzyme activator activity; enzyme binding; insulin receptor substrate binding; metal ion binding; non-membrane spanning protein tyrosine kinase activity; nucleotide binding; protein C-terminus bi