



Anti-CREB3L4 polyclonal antibody (DPABT-H20888)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Anti-CREB3L4 Polyclonal Antibody
Target	CREB3L4
Immunogen	Peptide
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Antigen Affinity Purified
Conjugate	Unconjugated
Applications	ELISA, WB, IHC
Molecular Weight	395aa,43kd
Cellular Localization	membrane
Format	Lyophilized Powder
Size	100 µl, 200 µl
Buffer	PBS with 0.02% Sodium Azide, 0.5%BSA,PH7.2
Preservative	0.02% Sodium Azide
Storage	-20 °C, Avoid freeze / thaw cycles

GENE INFORMATION

Gene Name	CREB3L4 cAMP responsive element binding protein 3-like 4 [Homo sapiens]
Official Symbol	CREB3L4
Synonyms	CREB3L4; cAMP responsive element binding protein 3-like 4; cyclic AMP-responsive element-binding protein 3-like protein 4; AlbZIP; ATCE1; CREB3; CREB4; hJAL; CREB-4; tisp40; attaching to CRE-like 1; cAMP responsive element binding protein 1; cAMP-responsi
Entrez Gene ID	148327
Protein Refseq	NP_001242907
UniProt ID	Q8TEY5
Chromosome Location	1q21.2
Pathway	Amphetamine addiction, organism-specific biosystem; Amphetamine addiction, conserved biosystem; Cholinergic synapse, organism-specific biosystem; Cocaine addiction, organism-specific biosystem; Cocaine addiction, conserved biosystem; Dopaminergic synapse,
Function	protein dimerization activity; sequence-specific DNA binding; sequence-specific DNA binding transcription factor activity;