



Anti-GNRH1 polyclonal antibody (DPABT-H21055)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Anti-GNRH1 Polyclonal Antibody
Target	GNRH1
Immunogen	Peptide
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Antigen Affinity Purified
Conjugate	Unconjugated
Applications	ELISA, IHC
Molecular Weight	92aa
Cellular Localization	Secreted
Format	Lyophilized Powder
Size	100 µl, 200 µl
Buffer	PBS with 0.02% Sodium Azide, 2%BSA,PH7.2
Preservative	0.02% Sodium Azide
Storage	-20 °C, Avoid freeze / thaw cycles

GENE INFORMATION

Gene Name	GNRH1 gonadotropin-releasing hormone2(luteinizing-releasing hormone) [Homo sapiens]
Official Symbol	GNRH1
Synonyms	GNRH1; gonadotropin-releasing hormone2(luteinizing-releasing hormone); GNRH, gonadotropin releasing hormone2(leutinizing releasing hormone) , GRH, LHRH; progonadoliberin-1; luliberin I; progonadoliberin I; GnRH-associated peptide 1; prolactin release-inh
Entrez Gene ID	2796
Protein Refseq	NP_000816
UniProt ID	P01148
Chromosome Location	8p21-p11.2
Pathway	Class A/1 (Rhodopsin-like receptors), organism-specific biosystem; G alpha (q) signalling events, organism-specific biosystem; GPCR downstream signaling, organism-specific biosystem; GPCR ligand binding, organism-specific biosystem; GnRH signaling pathway
Function	gonadotropin hormone-releasing hormone activity; hormone activity;