



Anti-FABP2 polyclonal antibody (DPABT-H21059)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Anti-FABP2 Polyclonal Antibody
Target	FABP2
Immunogen	Peptide
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Antigen Affinity Purified
Conjugate	Unconjugated
Applications	ELISA, IHC
Molecular Weight	132aa
Cellular Localization	Cytoplasmic
Format	Lyophilized Powder
Size	100 µl, 200 µl
Buffer	PBS with 0.02% Sodium Azide, 2%BSA,PH7.2
Preservative	0.02% Sodium Azide
Storage	-20 °C, Avoid freeze / thaw cycles

GENE INFORMATION

Gene Name	FABP2 fatty acid binding protein 2, intestinal [Homo sapiens]
Official Symbol	FABP2
Synonyms	FABP2; fatty acid binding protein 2, intestinal; fatty acid-binding protein, intestinal; I FABP; fatty acid-binding protein 2; intestinal-type fatty acid-binding protein; FABPI; I-FABP; MGC133132;
Entrez Gene ID	2169
Protein Refseq	NP_000125
UniProt ID	P12104
Chromosome Location	4q28-q31
Pathway	Fat digestion and absorption, organism-specific biosystem; Fat digestion and absorption, conserved biosystem; PPAR signaling pathway, organism-specific biosystem; PPAR signaling pathway, conserved biosystem;
Function	fatty acid binding; long-chain fatty acid transporter activity; transporter activity;