



Anti-KLK14 polyclonal antibody (DPABT-H21061)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Anti-KLK14 Polyclonal Antibody
Target	KLK14
Immunogen	Peptide
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Antigen Affinity Purified
Conjugate	Unconjugated
Applications	ELISA, IHC
Molecular Weight	251aa
Cellular Localization	Secreted; extracellular space.
Format	Lyophilized Powder
Size	100 µl, 200 µl
Buffer	PBS with 0.02% Sodium Azide, 2%BSA,PH7.2
Preservative	0.02% Sodium Azide
Storage	-20 °C, Avoid freeze / thaw cycles

GENE INFORMATION

Gene Name	KLK14 kallikrein-related peptidase 14 [Homo sapiens]
Official Symbol	KLK14
Synonyms	KLK14; kallikrein-related peptidase 14; kallikrein 14; kallikrein-14; KLK L6; hK14; kallikrein-like protein 6; KLK-L6;
Entrez Gene ID	43847
Protein Refseq	NP_071329
UniProt ID	A0A024R4M6
Chromosome Location	19q13.3-q13.4
Function	peptidase activity; serine-type endopeptidase activity; serine-type endopeptidase activity;