



Anti-PPP1R1B polyclonal antibody (DPABT-H21065)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Anti-PPP1R1B Polyclonal Antibody
Target	PPP1R1B
Immunogen	Peptide
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Antigen Affinity Purified
Conjugate	Unconjugated
Applications	ELISA, WB, IHC
Molecular Weight	204aa,32KD
Cellular Localization	Cytoplasm.
Format	Lyophilized Powder
Size	100 µl, 200 µl
Buffer	PBS with 0.02% Sodium Azide, 2%BSA,PH7.2
Preservative	0.02% Sodium Azide
Storage	-20 °C, Avoid freeze / thaw cycles

GENE INFORMATION

Gene Name	PPP1R1B protein phosphatase 1, regulatory (inhibitor) subunit 1B (dopamine and cAMP regulated phosphoprotein, DARPP-32) [Homo sapiens]
Official Symbol	PPP1R1B
Synonyms	PPP1R1B; protein phosphatase 1, regulatory (inhibitor) subunit 1B (dopamine and cAMP regulated phosphoprotein, DARPP-32); DARPP 32; dopamine and cAMP regulated phosphoprotein; FLJ20940; DARPP-32;
Entrez Gene ID	5503
Chromosome Location	17q12