



Anti-LRMP polyclonal antibody (DPABT-H21069)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Anti-LRMP Polyclonal Antibody
Target	LRMP
Immunogen	Peptide
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Antigen Affinity Purified
Conjugate	Unconjugated
Applications	ELISA, WB, IHC
Molecular Weight	499aa,
Cellular Localization	Processed lymphoid-restricted membrane protein: Cytoplasm By similarity.
Format	Lyophilized Powder
Size	100 µl, 200 µl
Buffer	PBS with 0.02% Sodium Azide, 2%BSA,PH7.2
Preservative	0.02% Sodium Azide
Storage	-20 °C, Avoid freeze / thaw cycles

GENE INFORMATION

Gene Name	LRMP lymphoid-restricted membrane protein [Homo sapiens]
Official Symbol	LRMP
Synonyms	LRMP; lymphoid-restricted membrane protein; JAW1; protein Jaw1;
Entrez Gene ID	4033
Protein Refseq	NP_001191055
UniProt ID	Q12912
Chromosome Location	12p12.1