



Anti-ARHGDIB polyclonal antibody (DPABT-H21070)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Anti-ARHGDIB Polyclonal Antibody
Target	ARHGDIB
Immunogen	Peptide
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Antigen Affinity Purified
Conjugate	Unconjugated
Applications	ELISA, WB, IHC
Molecular Weight	201aa,28KD
Cellular Localization	Cytoplasm.
Format	Lyophilized Powder
Size	100 µl, 200 µl
Buffer	PBS with 0.02% Sodium Azide, 2%BSA,PH7.2
Preservative	0.02% Sodium Azide
Storage	-20 °C, Avoid freeze / thaw cycles

GENE INFORMATION

Gene Name	ARHGDIB Rho GDP dissociation inhibitor (GDI) beta [Homo sapiens]
Official Symbol	ARHGDIB
Synonyms	ARHGDIB; Rho GDP dissociation inhibitor (GDI) beta; GDIA2, GDID4, RAP1GN1; rho GDP-dissociation inhibitor 2; Ly GDI; RhoGDI2; Rho GDI 2; rho-GDI beta; D4; GDIA2; GDID4; LYGDI; Ly-GDI; RAP1GN1;
Entrez Gene ID	397
Protein Refseq	NP_001166
UniProt ID	A0A024RAS5
Chromosome Location	12p12.3
Pathway	Caspase cascade in apoptosis, organism-specific biosystem; FAS pathway and Stress induction of HSP regulation, organism-specific biosystem; G13 Signaling Pathway, organism-specific biosystem; Neurotrophin signaling pathway, organism-specific biosystem; Ne
Function	GTPase activator activity; Rho GDP-dissociation inhibitor activity;