



Anti-ATF5 polyclonal antibody (DPABT-H21079)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Anti-ATF5 Polyclonal Antibody
Target	ATF5
Immunogen	Peptide
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Antigen Affinity Purified
Conjugate	Unconjugated
Applications	ELISA, IHC
Molecular Weight	282aa
Cellular Localization	Cytoplasm. Nucleus
Format	Lyophilized Powder
Size	100 µl, 200 µl
Buffer	PBS with 0.02% Sodium Azide, 2%BSA,PH7.2
Preservative	0.02% Sodium Azide
Storage	-20 °C, Avoid freeze / thaw cycles

GENE INFORMATION

Gene Name	ATF5 activating transcription factor 5 [Homo sapiens]
Official Symbol	ATF5
Synonyms	ATF5; activating transcription factor 5; cyclic AMP-dependent transcription factor ATF-5; transcription factor ATFx; cAMP-dependent transcription factor ATF-5; ATFX; HMFN0395; FLJ34666;
Entrez Gene ID	22809
Protein Refseq	NP_001180575
UniProt ID	A0A024QZG3
Chromosome Location	19q13.3
Pathway	Myometrial Relaxation and Contraction Pathways, organism-specific biosystem; Signaling events mediated by PRL, organism-specific biosystem;
Function	protein binding; protein dimerization activity; sequence-specific DNA binding; sequence-specific DNA binding transcription factor activity; transcription corepressor activity;