



Anti-DNM1 polyclonal antibody (DPABT-H21094)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Anti-DNM1 Polyclonal Antibody
Target	DNM1
Immunogen	Peptide
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Antigen Affinity Purified
Conjugate	Unconjugated
Applications	ELISA, IHC
Molecular Weight	864aa,97kd
Cellular Localization	Cytoplasm. Cytoplasm, cytoskeleton
Format	Lyophilized Powder
Size	100 µl, 200 µl
Buffer	PBS with 0.02% Sodium Azide, 2%BSA,PH7.2
Preservative	0.02% Sodium Azide
Storage	-20 °C, Avoid freeze / thaw cycles

GENE INFORMATION

Gene Name	DNM1 dynamin2[Homo sapiens]
Official Symbol	DNM1
Synonyms	DNM1; dynamin 1; DNM; dynamin-1;
Entrez Gene ID	1759
Protein Refseq	NP_001005336
UniProt ID	Q05193
Chromosome Location	9q34
Pathway	Axon guidance, organism-specific biosystem; Bacterial invasion of epithelial cells, organism-specific biosystem; Bacterial invasion of epithelial cells, conserved biosystem; CXCR3-mediated signaling events, organism-specific biosystem; CXCR4-mediated sign
Function	GTP binding; GTPase activity; hydrolase activity; identical protein binding; nucleotide binding; protein binding;