



Anti-eIF2C4 polyclonal antibody (DPABT-H21096)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Anti-EIF2C4 Polyclonal Antibody
Target	EIF2C4
Immunogen	Peptide
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Antigen Affinity Purified
Conjugate	Unconjugated
Applications	ELISA, IHC
Molecular Weight	861aa,97kd
Format	Lyophilized Powder
Size	100 µl, 200 µl
Buffer	PBS with 0.02% Sodium Azide, 2%BSA,PH7.2
Preservative	0.02% Sodium Azide
Storage	-20 °C, Avoid freeze / thaw cycles

GENE INFORMATION

Gene Name	EIF2C4 eukaryotic translation initiation factor 2C, 4 [Homo sapiens]
Official Symbol	EIF2C4
Synonyms	EIF2C4; eukaryotic translation initiation factor 2C, 4; protein argonaute-4; AGO4; argonaute 4; FLJ20033; KIAA1567; hAgo4; eIF2C 4; eIF-2C 4; argonaute4;
Entrez Gene ID	192670
Protein Refseq	NP_060099
UniProt ID	Q9HCK5
Chromosome Location	1p34
Pathway	Gene Expression, organism-specific biosystem; MicroRNA (miRNA) Biogenesis, organism-specific biosystem; Post-transcriptional Silencing By Small RNAs, organism-specific biosystem; Pre-NOTCH Expression and Processing, organism-specific biosystem; Pre-NOTCH
Function	RNA binding; protein binding;