



Anti-AQP10 polyclonal antibody (DPABT-H21099)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Anti-AQP10 Polyclonal Antibody
Target	AQP10
Immunogen	Peptide
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Antigen Affinity Purified
Conjugate	Unconjugated
Applications	ELISA, IHC
Molecular Weight	301aa,32KD
Cellular Localization	Membrane; Multi-pass membrane protein
Format	Lyophilized Powder
Size	100 µl, 200 µl
Buffer	PBS with 0.02% Sodium Azide, 2%BSA,PH7.2
Preservative	0.02% Sodium Azide
Storage	-20 °C, Avoid freeze / thaw cycles

GENE INFORMATION

Gene Name	AQP10 aquaporin 10 [Homo sapiens]
Official Symbol	AQP10
Synonyms	AQP10; aquaporin 10; aquaporin-10; AQP-10; aquaglyceroporin-10; small intestine aquaporin; AQPA_HUMAN;
Entrez Gene ID	89872
Protein Refseq	NP_536354
UniProt ID	Q96PS8
Chromosome Location	1q21.3
Pathway	Aquaporin-mediated transport, organism-specific biosystem; Passive Transport by Aquaporins, organism-specific biosystem; Transmembrane transport of small molecules, organism-specific biosystem.
Function	transporter activity;