



Anti-TRAK2 polyclonal antibody (DPABT-H21102)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Anti-TRAK2 Polyclonal Antibody
Target	TRAK2
Immunogen	Peptide
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Antigen Affinity Purified
Conjugate	Unconjugated
Applications	ELISA, IHC
Molecular Weight	914aa,101kDa
Format	Lyophilized Powder
Size	100 µl, 200 µl
Buffer	PBS with 0.02% Sodium Azide, 2%BSA,PH7.2
Preservative	0.02% Sodium Azide
Storage	-20 °C, Avoid freeze / thaw cycles

GENE INFORMATION

Gene Name	TRAK2 trafficking protein, kinesin binding 2 [Homo sapiens]
Official Symbol	TRAK2
Synonyms	TRAK2; trafficking protein, kinesin binding 2; ALS2CR3, amyotrophic lateral sclerosis 2 (juvenile) chromosome region, candidate 3; trafficking kinesin-binding protein 2; CALS C; gamma aminobutyric acid(A) receptor interacting factor; GRIF 1; KIAA0549; MI
Entrez Gene ID	66008
Protein Refseq	NP_055864
UniProt ID	O60296
Chromosome Location	2q33.1
Pathway	GABAergic synapse, organism-specific biosystem; GABAergic synapse, conserved biosystem; Metabolic pathways, organism-specific biosystem.
Function	GABA receptor binding; GABA receptor binding; enzyme binding; protein binding; receptor binding;