



Anti-ANLN polyclonal antibody (DPABT-H21109)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Anti-ANLN Polyclonal Antibody
Target	ANLN
Immunogen	Peptide
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Antigen Affinity Purified
Conjugate	Unconjugated
Applications	ELISA, IHC
Molecular Weight	1124 aa,124KD
Cellular Localization	Nucleus. Cytoplasm › cytoskeleton. Cytoplasm › cell cortex
Format	Lyophilized Powder
Size	100 µl, 200 µl
Buffer	PBS with 0.02% Sodium Azide, 2%BSA,PH7.2
Preservative	0.02% Sodium Azide
Storage	-20 °C, Avoid freeze / thaw cycles

GENE INFORMATION

Gene Name	ANLN anillin, actin binding protein [Homo sapiens]
Official Symbol	ANLN
Synonyms	ANLN; anillin, actin binding protein; anillin (Drosophila Scraps homolog), actin binding protein , anillin, actin binding protein (scraps homolog, Drosophila); actin-binding protein anillin; ANILLIN; scra; Scraps; DKFZp779A055;
Entrez Gene ID	54443
Protein Refseq	NP_061155
UniProt ID	Q9NQW6
Chromosome Location	7p15-p14
Function	actin binding;