



Anti-ADPRHL2 polyclonal antibody (DPABT-H21110)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Anti-ADPRHL2 Polyclonal Antibody
Target	ADPRHL2
Immunogen	Peptide
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Antigen Affinity Purified
Conjugate	Unconjugated
Applications	ELISA, IHC
Molecular Weight	363aa,39kd
Format	Lyophilized Powder
Size	200 µl
Buffer	PBS with 0.02% Sodium Azide, 2%BSA,PH7.2
Preservative	0.02% Sodium Azide
Storage	-20 °C, Avoid freeze / thaw cycles

GENE INFORMATION

Gene Name	ADPRHL2 ADP-ribosylhydrolase like 2 [Homo sapiens]
Official Symbol	ADPRHL2
Synonyms	ADPRHL2; ADP-ribosylhydrolase like 2; poly(ADP-ribose) glycohydrolase ARH3; ARH3; FLJ20446; ADP-ribosylhydrolase 3; protein ADP-ribosylarginine hydrolase-like protein 2; [Protein ADP-ribosylarginine] hydrolase-like protein 2;
Entrez Gene ID	54936
Protein Refseq	NP_060295
UniProt ID	Q9NX46
Chromosome Location	1p34.3
Function	hydrolase activity; metal ion binding; poly(ADP-ribose) glycohydrolase activity;