



Anti-ADAMTS7 polyclonal antibody (DPABT-H21115)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Anti-ADAMTS7 Polyclonal Antibody
Target	ADAMTS7
Immunogen	Peptide
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Serum
Conjugate	Unconjugated
Applications	ELISA, WB
Molecular Weight	997aa,110kd
Cellular Localization	Secreted, extracellular space, extracellular matrix (By similarity)
Format	Lyophilized Powder
Size	200 µl
Buffer	Serum with 0.02% Sodium Azide
Preservative	0.02% Sodium Azide
Storage	-20 °C, Avoid freeze / thaw cycles

GENE INFORMATION

Gene Name	ADAMTS7 ADAM metallopeptidase with thrombospondin type2motif, 7 [Homo sapiens]
Official Symbol	ADAMTS7
Synonyms	ADAMTS7; ADAM metallopeptidase with thrombospondin type2motif, 7; a disintegrin like and metalloprotease (reprolysin type) with thrombospondin type2motif, 7; A disintegrin and metalloproteinase with thrombospondin motifs 7; a disintegrin and metalloprotea
Entrez Gene ID	11173
Protein Refseq	NP_055087
UniProt ID	Q9UFZ4
Chromosome Location	15pter-qter
Function	metal ion binding; metalloendopeptidase activity; metallopeptidase activity; peptidase activity; zinc ion binding;