



Anti-MST1R polyclonal antibody (DPABT-H21122)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Anti-MST1R Polyclonal Antibody
Target	MST1R
Immunogen	Peptide
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Antigen Affinity Purified
Conjugate	Unconjugated
Applications	ELISA, IHC
Molecular Weight	1400aa,152kd
Cellular Localization	Single-pass type I membrane protein.
Format	Lyophilized Powder
Size	200 µl
Buffer	PBS with 0.02% Sodium Azide, 2%BSA,PH7.2
Preservative	0.02% Sodium Azide
Storage	-20 °C, Avoid freeze / thaw cycles

GENE INFORMATION

Gene Name	MST1R macrophage stimulating 1 receptor (c-met-related tyrosine kinase) [Homo sapiens]
Official Symbol	MST1R
Synonyms	MST1R; macrophage stimulating 1 receptor (c-met-related tyrosine kinase); PTK8, PTK8 protein tyrosine kinase 8 , RON; macrophage-stimulating protein receptor; CD136; CDw136; p185-Ron; MSP receptor; RON variant E2E3; MST1R variant RON30; MST1R variant RON62; soluble RON variant 1; soluble RON variant 2; soluble RON variant 3; soluble RON variant 4; c-met-related tyrosine kinase; PTK8 protein tyrosine kinase 8; RON; PTK8;
Entrez Gene ID	4486
Protein Refseq	NP_001231866
UniProt ID	B4E058
Chromosome Location	3p21
Pathway	Alpha6-Beta4 Integrin Signaling Pathway, organism-specific biosystem; a6b1 and a6b4 Integrin signaling, organism-specific biosystem; amb2 Integrin signaling, organism-specific biosystem;
Function	ATP binding; enzyme binding; macrophage colony-stimulating factor receptor activity; nucleotide binding; protein binding; receptor activity;