



Anti-TNFRSF19 polyclonal antibody (DPABT-H21146)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Anti-TNFRSF19 Polyclonal Antibody
Target	TNFRSF19
Immunogen	Peptide
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Serum
Conjugate	Unconjugated
Applications	ELISA, WB
Molecular Weight	423aa,46kd
Cellular Localization	Membrane; Single-pass type I membrane protein (Probable)
Format	Lyophilized Powder
Size	200 µl
Buffer	Serum with 0.02% Sodium Azide
Preservative	0.02% Sodium Azide
Storage	-20 °C, Avoid freeze / thaw cycles

GENE INFORMATION

Gene Name	TNFRSF19 tumor necrosis factor receptor superfamily, member 19 [Homo sapiens]
Official Symbol	TNFRSF19
Synonyms	TNFRSF19; tumor necrosis factor receptor superfamily, member 19; tumor necrosis factor receptor superfamily member 19; TAJ; TAJ alpha; toxicity and JNK inducer; TRADE; TROY; TAJ-alpha;
Entrez Gene ID	55504
Protein Refseq	NP_001191387
UniProt ID	Q9NS68
Chromosome Location	13q12.11-q12.3
Pathway	Cytokine-cytokine receptor interaction, organism-specific biosystem; Cytokine-cytokine receptor interaction, conserved biosystem;
Function	binding; receptor activity; tumor necrosis factor-activated receptor activity;