



Rabbit anti-RELA Polyclonal antibody (DPABT-H21147)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Anti-RELA Polyclonal Antibody
Target	RELA
Immunogen	Peptide
Isotype	IgG
Source/Host	Rabbit
Purification	Serum
Conjugate	Unconjugated
Applications	ELISA, WB, IHC
Molecular Weight	551aa,65kd
Cellular Localization	Nucleus. Cytoplasm
Format	Lyophilized Powder
Size	200 µl
Buffer	Serum with 0.02% Sodium Azide
Preservative	0.02% Sodium Azide
Storage	-20 °C, Avoid freeze / thaw cycles

GENE INFORMATION

Gene Name	RELA v-rel reticuloendotheliosis viral oncogene homolog A (avian) [Homo sapiens]
Official Symbol	RELA
Synonyms	RELA; v-rel reticuloendotheliosis viral oncogene homolog A (avian); NFKB3, nuclear factor of kappa light polypeptide gene enhancer in B cells 3; transcription factor p65; p65; NF-kappa-B p65delta3; nuclear factor NF-kappa-B p65 subunit; nuclear factor of kappa light polypeptide gene enhancer in B-cells 3; NFKB3; MGC131774;
Entrez Gene ID	5970
Protein Refseq	NP_001138610
UniProt ID	Q04206
Chromosome Location	11q13
Pathway	Activated TLR4 signalling, organism-specific biosystem; Activation of NF-kappaB in B Cells, organism-specific biosystem; Acute myeloid leukemia, organism-specific biosystem; Acute myeloid leukemia, conserved biosystem; Adaptive Immune System, organism-specific biosystem; Adipocytokine signaling pathway, organism-specific biosystem; Adipocytokine signaling pathway, conserved biosystem;
Function	DNA binding; NF-kappaB binding; activating transcription factor binding; ankyrin repeat binding; chromatin binding; identical protein binding; identical protein binding; phosphate ion binding; protein N-terminus binding; protein binding; protein kinase bi