



Anti-EREG polyclonal antibody (DPABT-H21153)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Anti-EREG Polyclonal Antibody
Target	EREG
Immunogen	Peptide
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Antigen Affinity Purified
Conjugate	Unconjugated
Applications	ELISA, WB, IHC
Molecular Weight	169aa,19kd
Cellular Localization	Secreted
Format	Lyophilized Powder
Size	100 µl, 200 µl
Buffer	PBS with 0.02% Sodium Azide, 0.5%BSA,PH7.2
Preservative	0.02% Sodium Azide
Storage	-20 °C, Avoid freeze / thaw cycles

GENE INFORMATION

Gene Name	EREG epiregulin [Homo sapiens]
Official Symbol	EREG
Synonyms	EREG; epiregulin; proepiregulin; ER;
Entrez Gene ID	2069
Protein Refseq	NP_001423
UniProt ID	O14944
Chromosome Location	4q21.21
Pathway	ErbB receptor signaling network, organism-specific biosystem; ErbB signaling pathway, organism-specific biosystem; ErbB signaling pathway, organism-specific biosystem; ErbB signaling pathway, conserved biosystem; ErbB4 signaling events, organism-specific biosystem; GRB2 events in ERBB2 signaling, organism-specific biosystem; Nuclear signaling by ERBB4, organism-specific biosystem;
Function	epidermal growth factor receptor binding; epidermal growth factor receptor binding; growth factor activity; protein binding;