



Anti-FASLG polyclonal antibody (DPABT-H21161)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Anti-FASLG Polyclonal Antibody
Target	FASLG
Immunogen	Peptide
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Serum
Conjugate	Unconjugated
Applications	ELISA, WB, IHC
Molecular Weight	281aa,31kd
Cellular Localization	Cell membrane; Single-pass type II membrane protein. Secreted
Format	Lyophilized Powder
Size	200 µl
Buffer	Serum with 0.02% Sodium Azide
Preservative	0.02% Sodium Azide
Storage	-20 °C, Avoid freeze / thaw cycles

GENE INFORMATION

Gene Name	FASLG Fas ligand (TNF superfamily, member 6) [Homo sapiens]
Official Symbol	FASLG
Synonyms	FASLG; Fas ligand (TNF superfamily, member 6); APT1LG1, TNFSF6, tumor necrosis factor (ligand) superfamily, member 6; tumor necrosis factor ligand superfamily member 6; CD178; FasL; APTL; CD95 ligand; fas antigen ligand; apoptosis antigen ligand; apoptos
Entrez Gene ID	356
Protein Refseq	NP_000630
UniProt ID	P48023
Chromosome Location	1q23
Pathway	Activation of Pro-Caspase 8, organism-specific biosystem; African trypanosomiasis, organism-specific biosystem; African trypanosomiasis, conserved biosystem; Allograft rejection, organism-specific biosystem; Allograft rejection, conserved biosystem; Apopt
Function	cytokine activity; protein binding; receptor binding; tumor necrosis factor receptor binding;