



Anti-TNK2 polyclonal antibody (DPABT-H21165)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Anti-TNK2 Polyclonal Antibody
Target	TNK2
Immunogen	Peptide
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse
Purification	Serum
Conjugate	Unconjugated
Applications	ELISA, WB
Molecular Weight	1055aa,117kd
Cellular Localization	Cell membrane
Format	Lyophilized Powder
Size	200 µl
Buffer	Serum with 0.02% Sodium Azide
Preservative	0.02% Sodium Azide
Storage	-20 °C, Avoid freeze / thaw cycles

GENE INFORMATION

Gene Name	TNK2 tyrosine kinase, non-receptor, 2 [Homo sapiens]
Official Symbol	TNK2
Synonyms	TNK2; tyrosine kinase, non-receptor, 2; activated CDC42 kinase 1; ACK; ACK1; activated Cdc42 associated kinase 1; p21cdc42Hs; activated p21cdc42Hs kinase; activated Cdc42-associated kinase 1; tyrosine kinase non-receptor protein 2; ACK-1; FLJ44758; FLJ455
Entrez Gene ID	10188
Protein Refseq	NP_001010938
UniProt ID	Q07912
Chromosome Location	3q29
Pathway	CDC42 signaling events, organism-specific biosystem; EGFR1 Signaling Pathway, organism-specific biosystem; Focal Adhesion, organism-specific biosystem; G13 Signaling Pathway, organism-specific biosystem;
Function	ATP binding; GTPase inhibitor activity; metal ion binding; non-membrane spanning protein tyrosine kinase activity; nucleotide binding; protein binding; protein serine/threonine kinase activity; protein serine/threonine/tyrosine kinase activity; protein ty