



Anti-ANXA10 polyclonal antibody (DPABT-H21170)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Anti-ANXA10 Polyclonal Antibody
Target	ANXA10
Immunogen	Peptide
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Antigen Affinity Purified
Conjugate	Unconjugated
Applications	ELISA, WB
Molecular Weight	324aa,37kd
Format	Lyophilized Powder
Size	200 µl
Buffer	PBS with 0.02% Sodium Azide, 2%BSA,PH7.2
Preservative	0.02% Sodium Azide
Storage	-20 °C, Avoid freeze / thaw cycles

GENE INFORMATION

Gene Name	ANXA10 annexinB20[Homo sapiens]
Official Symbol	ANXA10
Synonyms	ANXA10; annexin A10; ANX14; annexin 14; annexin-10; annexin-14;
Entrez Gene ID	11199
Protein Refseq	NP_009124
UniProt ID	Q9UJ72
Chromosome Location	4q32.3
Function	calcium ion binding; calcium-dependent phospholipid binding;