



Rabbit anti-CapG Polyclonal antibody (DPABT-H21172)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Anti-CapG Polyclonal Antibody
Target	CapG
Immunogen	Peptide
Isotype	IgG
Source/Host	Rabbit
Purification	Antigen Affinity Purified
Conjugate	Unconjugated
Applications	ELISA, WB, IHC
Molecular Weight	348aa,39kd
Cellular Localization	Cytoplasm. Nucleus. Melanosome
Format	Lyophilized Powder
Size	200 µl
Buffer	PBS with 0.02% Sodium Azide, 2%BSA,PH7.2
Preservative	0.02% Sodium Azide
Storage	-20 °C, Avoid freeze / thaw cycles

GENE INFORMATION

Gene Name	CAPG capping protein (actin filament), gelsolin-like [Homo sapiens]
Official Symbol	CapG
Synonyms	CAPG; capping protein (actin filament), gelsolin-like; AFCP; macrophage-capping protein; macrophage capping protein; MCP; gelsolin-like capping protein; actin regulatory protein CAP-G; actin-regulatory protein CAP-G;
Entrez Gene ID	822
Protein Refseq	NP_001243068
UniProt ID	B2R9S4
Chromosome Location	2p11.2
Function	actin binding;