



Rabbit anti-ATXN7 Polyclonal antibody (DPABT-H21174)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Anti-ATXN7 Polyclonal Antibody
Target	ATXN7
Immunogen	Peptide
Isotype	IgG
Source/Host	Rabbit
Purification	Purified
Conjugate	Unconjugated
Applications	ELISA, WB
Molecular Weight	892aa,95kd
Cellular Localization	Nucleus
Format	Lyophilized Powder
Size	200 µl
Buffer	PBS with 0.02% Sodium Azide, PH7.2
Preservative	0.02% Sodium Azide
Storage	-20 °C, Avoid freeze / thaw cycles

GENE INFORMATION

Gene Name	ATXN7 ataxin 7 [Homo sapiens]
Official Symbol	ATXN7
Synonyms	ATXN7; ataxin 7; SCA7, spinocerebellar ataxia 7 (olivopontocerebellar atrophy with retinal degeneration); ataxin-7; ADCAI1; OPCA3; spinocerebellar ataxia type 7 protein; SCA7; FLJ117787;
Entrez Gene ID	6314
Protein Refseq	NP_000324
UniProt ID	O15265
Chromosome Location	3p21.1-p12
Function	chromatin binding; protein binding; zinc ion binding;