



Anti-TNIP1 polyclonal antibody (DPABT-H21176)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Anti-TNIP1 Polyclonal Antibody
Target	TNIP1
Immunogen	Peptide
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Purified
Conjugate	Unconjugated
Applications	ELISA, WB
Molecular Weight	441aa,50kd
Cellular Localization	Nucleus
Format	Lyophilized Powder
Size	200 µl
Buffer	PBS with 0.02% Sodium Azide, PH7.2
Preservative	0.02% Sodium Azide
Storage	-20 °C, Avoid freeze / thaw cycles

GENE INFORMATION

Gene Name	TNIP1 TNFAIP3 interacting protein2[Homo sapiens]
Official Symbol	TNIP1
Synonyms	TNIP1; TNFAIP3 interacting protein 1; TNFAIP3-interacting protein 1; ABIN 1; KIAA0113; NAF1; Nef associated factor2SNP; VAN; virion associated nuclear shuttling protein; Nef-associated factor2SNP; HIV-1 Nef-interacting protein; virion-associated nuclear s
Entrez Gene ID	10318
Protein Refseq	NP_001239314
UniProt ID	B7Z8K2
Chromosome Location	5q32-q33.1
Pathway	EGFR1 Signaling Pathway, organism-specific biosystem; TNF-alpha/NF-kB Signaling Pathway, organism-specific biosystem;
Function	protein binding;