



Anti-AP3D1 polyclonal antibody (DPABT-H21178)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Anti-AP3D1 Polyclonal Antibody
Target	AP3D1
Immunogen	Peptide
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Serum
Conjugate	Unconjugated
Applications	ELISA, WB
Molecular Weight	1153aa,130kd
Cellular Localization	Golgi apparatus membrane; Peripheral membrane protein; Cytoplasmic side
Format	Lyophilized Powder
Size	200 µl
Buffer	Serum with 0.02% Sodium Azide
Preservative	0.02% Sodium Azide
Storage	-20 °C, Avoid freeze / thaw cycles

GENE INFORMATION

Gene Name	AP3D1 adaptor-related protein complex 3, delta2subunit [Homo sapiens]
Official Symbol	AP3D1
Synonyms	AP3D1; adaptor-related protein complex 3, delta2subunit; AP-3 complex subunit delta-1; ADTD; delta adaptin; adaptin, delta; AP-3 complex subunit delta; delta-adaptin, partial CDS; AP-3 complex delta subunit, partial CDS; adapter-related protein complex4su
Entrez Gene ID	8943
Protein Refseq	NP_001070991
UniProt ID	O14617
Chromosome Location	19p13.3
Pathway	Lysosome, organism-specific biosystem; Lysosome, conserved biosystem.
Function	binding; protein transporter activity; transporter activity;