



Rabbit anti-HLX Polyclonal antibody (DPABT-H21184)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Anti-HLX Polyclonal Antibody
Target	HLX
Immunogen	Peptide
Isotype	IgG
Source/Host	Rabbit
Purification	Serum
Conjugate	Unconjugated
Applications	ELISA, WB
Molecular Weight	488aa,51kd
Cellular Localization	Nucleus
Format	Lyophilized Powder
Size	200 µl
Buffer	Serum with 0.02% Sodium Azide
Preservative	0.02% Sodium Azide
Storage	-20 °C, Avoid freeze / thaw cycles

GENE INFORMATION

Gene Name	HLX H2.0-like homeobox [Homo sapiens]
Official Symbol	HLX
Synonyms	HLX; H2.0-like homeobox; H2.0 (Drosophila) like homeo box2, H2.0 like homeobox2, H2.0 like homeobox2(Drosophila) , HLX1; H2.0-like homeobox protein; HB24; H2.0-like homeobox 1; H2.0-like homeo box 1; H2.0-like homeo box-1; homeobox protein HB24; homeobo
Entrez Gene ID	3142
Protein Refseq	NP_068777
UniProt ID	Q14774
Chromosome Location	1q41
Pathway	IL12-mediated signaling events, organism-specific biosystem;
Function	sequence-specific DNA binding; transcription regulatory region sequence-specific DNA binding;